

Product Datasheet

SIAH1/2 Peptide (orb224377)

Description	The peptide is used to block SIAH1/2 antibody (orb224178) reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SIAH1; SIAH2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224178. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q920M9 , Q06986 , O43255 , Q8R4T2 , Q8IUQ4
Entrez	6478 , 140593 , 140941 , 6477 , 20439
Expiration Date	12 months from date of receipt.