

Product Datasheet

SP3/4 Peptide (orb224378)

Description	The peptide is used to block SP3/4 antibody (orb224179) reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SP3; SP4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224179. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02447 , O70494 , Q02446 , Q62445
Entrez	6670 , 6671 , 20687 , 20688
Expiration Date	12 months from date of receipt.