


Product Datasheet

CYP3A4 antibody (orb225118)

Description	Mouse monoclonal antibody to CYP3A4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	ELISA: 1-10000, WB: 1-500-1:2000, IHC-P: 1-200-1:1000, IF/ICC: 1-200-1:1000, FACS: 1-200-1:400
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant pro...