

Product Datasheet

SMCO3 antibody (Biotin) (orb22536)

Description	SMCO3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Single-pass membrane and coiled-coil domain-containing protein 3 protein (1-154AA)
Target	SMCO3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A2RU48
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.