



## Product Datasheet

### PDLIM1 antibody (orb225557)


<b>Description</b>	Rabbit polyclonal antibody to PDLIM1
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Human, Mouse, Rabbit, Rat, Sheep, Xenopus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IF, IHC-P, WB
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB: 1:1000-3000, IHC-P: 1:200, IF/ICC: 1:200-500, ELISA: 1:20000-40000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Jurkat cells, C...



Western blot analysis of extracts from r...



Western blot analysis of extracts from 2...