

Product Datasheet

G CSF antibody (Biotin) (orb22589)

Description	Mouse monoclonal antibody against G CSF conjugated to Biotin
Reactivity	Human
Conjugation	Unconjugated
Immunogen	r.Human G-CSF
Preservatives	1 mg/ml in PBS (after reconstitution).
Storage	Storage Procedures: In lyophilized form, for long periods, store at 4°C in a dry environment. After reconstitution, if not intended for use within a month, aliquot and store at -20°C
Note	For research use only
Application notes	Applications: Direct ELISA, Western Blot, Intra-cellular staining. Titer: For intra-cellular staining use 10µl per 1,000,000 cells
Isotype	Mouse IgG
Clonality	Monoclonal
Solubility (25°C)	Reconstitute with sterile H2O. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
Expiration Date	12 months from date of receipt.