

Product Datasheet

CD90 Thy-1.2 Antibody (orb22707)

Description	Mouse monoclonal antibody to CD90 Thy-1.2
Reactivity	Mouse
Conjugation	Unconjugated
Immunogen	Purified mouse LN T cells
Preservatives	1mg/ml in PBS (after reconstitution).
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: Cytotoxic and staining antibody. For staining, use 10µl/1,000,000 cells. Exact titer for cytotoxicity should be determined by the investigator
Isotype	Rat IgM
Clonality	Monoclonal
Clone Number	mThy-1.2
Solubility (25°C)	Reconstitute with of H2O. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)