

## Product Datasheet

### Dengue type 2 antibody (orb22709)

<b>Description</b>	Mouse monoclonal antibody to Dengue type 2
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Dengue-infected cell lysate
<b>Preservatives</b>	1mg/ml in PBS (after reconstitution).
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Note: Minimal cross-reactivity with other Dengue subtypes. Titer: By intracellular staining of Dengue-infected cells, 1:1000 dilution
<b>Isotype</b>	mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Solubility (25°C)</b>	Reconstitute with H2O. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
<b>Expiration Date</b>	12 months from date of receipt.