

Product Datasheet

HCV NS3 antibody (orb22713)

Description	Mouse monoclonal antibody to HCV NS3
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Recombinant HCV NS-3 (genotype 1a)
Preservatives	Lyophilized from a solution containing no additives.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Monoclonal Antibodies, Anti-Viral
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.