

Product Datasheet

Influenza-A H3N2 antibody (orb22721)

Description	Mouse monoclonal antibody to Influenza-A H3N2
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Influenza A hemagglutinin H3N2
Preservatives	In PBS, pH 7.4 and 0.1 % NaN ₃ .
Form/Appearance	Sterile Filtered clear solution.
Storage	Should be stored at 4C. Do not freeze.
Note	For research use only
Application notes	Applications: Influenza A haemagglutinin H3N2 immunodetection in direct or indirect ELISA, and Western Blotting
Isotype	mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.