

Product Datasheet

HIV-1 gp41 antibody (orb22731)

Description	Rabbit polyclonal antibody to HIV-1 gp41
Reactivity	Virus
Conjugation	Unconjugated
Form/Appearance	Sterile Filtered solution.
Storage	Storage Procedures: Store at -20°C. Stability / Shelf Life: One year frozen. One month in solution at 4°C
Note	For research use only
Application notes	Polyclonal Antibodies, Anti-Viral
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.