


Product Datasheet

His-Tag Rabbit pAb antibody (orb227747)

Description	Rabbit polyclonal antibody to His-Tag Rabbit pAb
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:3000-10000, ELISA 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of human cell line...