

Product Datasheet

RORyt inverse agonist 31 (orb2277919)

Description	RORyt inverse agonist 31 (14g) is a potent antagonist of the retinoic acid receptor-related orphan receptor γ t (ROR γ t), with an inhibitory concentration (IC ₅₀) of 0.428 μ M. It has shown efficacy in reducing Imiquimod-induced psoriasis severity in murine models [1].
Storage	-20°C
Note	For research use only
Formula	$C_{23}H_{15}Cl_2F_3N_4O_3S$
Purity	98.00%
MW	555.36
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)