

Product Datasheet

Tubulin polymerization-IN-45 (orb2278019)

Description	Tubulin polymerization-IN-45, a tubulin-targeting agent, inhibits tubulin polymerization by binding to tubulin's colchicine site, inducing apoptotic cell death in hepatocellular cancer (HCC) cells [1].
Storage	-20°C
Note	For research use only
Formula	$C_{20}H_{18}N_4O_3$
Purity	98.00%
MW	362.38
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)