

Product Datasheet

IHMT-PI3K-455 (orb2278369)

Description	IHMT-PI3K-455 (Compound 15u) is a potent, selective PI3K γ / δ dual inhibitor with IC50 values of 7.1 nM for PI3K γ and 0.57 nM for PI3K δ , exhibiting oral activity. This compound inhibits AKT phosphorylation and promotes tumor regression by recruiting and activating CD8+ T cells. It is utilized in cancer research [1].
Storage	-20°C
Note	For research use only
Formula	$C_{26}H_{21}F_2N_7O_3$
Purity	98.00%
MW	517.49
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)