

Product Datasheet

Epcoritamab (orb2278696)

Description	Epcoritamab (GEN3013) is a novel subcutaneous CD3xCD20 bispecific T-cell binding antibody that activates T-cells and directs them to kill malignant CD20(+) B-cells for the study of relapsed or refractory large B-cell lymphoma.
Conjugation	Unconjugated
Target	Epcoritamab
Concentration	6.46 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Purity	98.00%
CAS Number	2134641-34-0
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)