

---

## Product Datasheet

**Mouse ASGR1 Protein, Fc Tag (MALS & SPR  
verified) (orb2279543)**

**Description**

Mouse ASGR1 Protein, Fc Tag is

(HEK293). It contains AA Gln 61 - Asn ...

**Reactivity**

Mouse

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

ASGR1

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

N-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 52.2 kDa. The protein migrates as 38-43 kDa and 55-65 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**MW**

52.2 kDa

**Uniprot ID****P34927****Expiration Date**

6 months from date of receipt.