
Product Datasheet

Human parainfluenza viruses Pre-fusion glycoprotein F0, His Tag (MALS verified) (orb2279544)

Description

Human parainfluenza viruses Pre-fusion

	human 293 cells (H...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Pre-Fusion glycoprotein F0
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.8 kDa. The protein migrates as 58-65 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
MW	57.8 kDa
Uniprot ID	AGW51052.1 (Q162C,L168C,I213C,G230C,A463V,I474Y)
Expiration Date	6 months from date of receipt.