
Product Datasheet

Influenza B (B/Singapore/INFTT-16-0610/2016)

Hemagglutinin (HA) Protein, His Tag (orb2279545)

Description

Influenza B (B/Singapore/INFTT-16-

Protein, His Tag is expressed from h...

Reactivity

Virus

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Hemagglutinin (HA)

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 59.3 kDa. The protein migrates as 55 kDa and 85-110 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

MW

59.3 kDa

Uniprot ID

EPI1843276

Expiration Date

6 months from date of receipt.