

---

## Product Datasheet

**Biotinylated Mouse Transferrin R / CD71 Protein,  
,His Tag (MALS verified) (orb2279546)**

**Description****Biotinylated Mouse Transferrin**

expressed from human 293 cells (HEK293)...

**Reactivity**

Mouse

**Endotoxins**

1.0 EU per µg

**Conjugation**

Biotin

**Target**

Transferrin R

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

N-Avi & N-10×His

**Note**

For research use only

**Application notes**

This protein carries an Avi tag (A) at the N-terminus, followed by a polyhistidine tag. The protein has a calculated MW of 79.5 kDa. The protein migrates as 85-95 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**MW**

79.5 kDa

**Uniprot ID**

[Q62351](#)

**Expiration Date**

6 months from date of receipt.