

---

## Product Datasheet

**Biotinylated Rat LRRC15 / LIB Protein, His, (MALS verified) (orb2279548)**

**Description**

Biotinylated Rat LRRC15 Protein,

cells (HEK293). It contains AA Tyr...

**Reactivity**

Rat

**Endotoxins**

1.0 EU per µg

**Conjugation**

Biotin

**Target**

LRRC15

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His &amp; C-Avi

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 61 kDa. The protein migrates as 70-90 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**MW**

61 kDa

**Uniprot ID**
[Q8R5M3](#)
**Expiration Date**

6 months from date of receipt.