

---

## Product Datasheet

**Human IL-18BP Protein, His, (SPR verified)  
(orb2279552)**

**Description**

Human IL-18BP Protein, His, is

(HEK293). It contains AA Thr 31 - Gly  
...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Biotin

**Target**

IL-18BP

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in  
lyophilized state;. -70°C for 3 months  
under sterile conditions after  
reconstitution. For long term storage,  
the product should be stored at  
lyophilized state at -20°C or lower

**Tag**

C-10×His &amp; C-Avi

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag  
at the C-terminus, followed by an Avi  
tag (). The protein has a calculated  
MW of 21.2 kDa. The protein migrates  
as 33-55 kDa under reducing (R)  
condition (SDS-PAGE) due to  
glycosylation.

**Purity**

90%

**MW**

21.2 kDa

**Uniprot ID**

NP\_001138529.1

**Expiration Date**

6 months from date of receipt.