

---

## Product Datasheet

**Human ANGPTL7 / CDT6 Protein, His Tag  
(orb2279553)**

**Description**

Human ANGPTL7 Protein, His Tag is

(HEK293). It contains AA Gln 27 - P...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

ANGPTL7

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.0 kDa. The protein migrates as 48-55 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**MW**

39.0 kDa

**Uniprot ID****043827****Expiration Date**

6 months from date of receipt.