
Product Datasheet

Human VEGF-C / Flt4-L Protein, His Tag (MALS verified) (orb2279554)

Description

Human VEGF-C Protein, His Tag is

(HEK293). It contains AA Thr 103 - A...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

VEGF-C

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 31 kDa. The protein migrates as 40-45 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

MW

31 kDa

Uniprot ID

AAH35212.1

Expiration Date

6 months from date of receipt.