

---

## Product Datasheet

**Human IL-17A / CTLA-8 Protein, His Tag (MALS verified) (orb2279555)**

**Description**

Human IL-17A Protein, His Tag is

(HEK293). It contains AA Gly 24 - Al...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

IL-17A

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

N-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 17.0 kDa. The protein migrates as 17-20 kDa and 22-25 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**MW**

17.0 kDa

**Uniprot ID**[Q16552](#)**Expiration Date**

6 months from date of receipt.