

---

## Product Datasheet

**Cynomolgus CDCP1 Protein, His Tag (MALS verified) (orb2279557)**

**Description**

Cynomolgus CDCP1 Protein, His Tag

(HEK293). It contains AA Phe 30 ...

**Reactivity**

Monkey

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

CDCP1

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 73.9 kDa. The protein migrates as 95-120 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**MW**

73.9 kDa

**Uniprot ID**

XP\_005546930.2

**Expiration Date**

6 months from date of receipt.