
Product Datasheet

**Mouse CDCP1 Protein, His Tag (MALS verified)
(orb2279558)**

Description

Mouse CDCP1 Protein, His Tag is

(HEK293). It contains AA Ser 30 - Leu...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CDCP1

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 73.7 kDa. The protein migrates as 90-95 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

MW

73.7 kDa

Uniprot ID
[Q5U462](#)
Expiration Date

6 months from date of receipt.