

## Product Datasheet

### Human IL-8 Protein (orb2279561)

<b>Description</b>	Human IL-8 Protein, premium grade is expressed from human 293 cells (HEK293). It contains AA Ser 28 - Ser 99 (Accession # P10145-1).
<b>Reactivity</b>	Human
<b>Endotoxins</b>	0.01 EU per $\mu$ g
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IL-8, CXCL8
<b>Preservatives</b>	Powder
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C to -70°C for 24 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
<b>Tag</b>	Native
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries no "tag". The protein has a calculated MW of 8.4 kDa. The protein migrates as 9 kDa $\pm$ 3 kDa under reducing (R) condition, and 11 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.
<b>Purity</b>	95%
<b>MW</b>	8.4 kDa
<b>Uniprot ID</b>	<b>P10145</b>
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)