
Product Datasheet

**Human IL-8 / CXCL8 Protein, premium grade
(orb2279561)**

Description

Human IL-8 Protein, premium grade

(HEK293). It contains AA Ser 28 ...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

IL-8

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 24 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 8.4 kDa. The protein migrates as 9 kDa±3 kDa under reducing (R) condition, and 11 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

Purity

95%

MW

8.4 kDa

Uniprot ID
P10145
Expiration Date

6 months from date of receipt.