
Product Datasheet

Human Thrombopoietin / TPO Protein, premium grade (orb2279562)

Description

Human Thrombopoietin Protein,

human 293 cells (HEK293). It contains ...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

Thrombopoietin

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 24 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 35.5 kDa. The protein migrates as 75 kDa±5 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

MW

35.5 kDa

Uniprot ID**P40225****Expiration Date**

6 months from date of receipt.