

Product Datasheet

Human IL-9 Protein, premium grade (orb2279564)

Description	Human IL-9 Protein, premium grade is expressed from human 293 cells (HEK293). It contains AA Gln 19 ...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	IL-9
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 24 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 14.1 kDa. The protein migrates as 30 kDa±3 kDa under reducing (R) condition, and 30 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.
Purity	95%
MW	14.1 kDa
Uniprot ID	P15248
Expiration Date	6 months from date of receipt.