

## Product Datasheet

### L343 (orb2279868)

<b>Description</b>	L-343 is an ionizable cationic lipidoid and can be used to synthetic liposomes for systemic delivery...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{51}H_{95}NO_6$
<b>MW</b>	818.3
<b>CAS Number</b>	[1443520-22-6]
<b>SMILES</b>	<chem>O=C(OC(C)(C)CCCCCCC/C=C\CCCCCCCC(OC(CCCN(C)C)=O)CCCCC CC/C=C\CCCCCCCC(OC(C)(C)C)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.