

## Product Datasheet

# Recombinant Monkeypox Virus (MPXV) A29L Monoclonal Antibody (orb2280050)

<b>Description</b>	Monkeypox virus (MPXV), which can cause monkeypox (MPX), is a zoonotic enveloped dsDNA virus that be...
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	LF
<b>Form/Appearance</b>	Purified,Liquid
<b>Storage</b>	For long term storage, the product should be stored at -20°C or lower.
<b>Note</b>	For research use only
<b>Application notes</b>	Recommended Applications: LFIA
<b>Sample Types</b>	Test sample type: Rash acne blister exudate / Nasopharyngeal swabs / serum. Test sample source: Human
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95% pure (SDS-PAGE).
<b>MW</b>	150 kDa
<b>Expiration Date</b>	12 months from date of receipt.