

Product Datasheet

Recombinant Monkeypox Virus (MPXV) M1R Monoclonal Antibody (orb2280053)

Description	Monkeypox virus (MPXV), which can cause monkeypox (MPX), is a zoonotic enveloped dsDNA virus that be...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA
Sample Types	Test sample type: Rash acne blister exudate / Nasopharyngeal swabs / serum. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.