


Product Datasheet

GRP78 antibody (orb228011)


Description	Rabbit polyclonal antibody to GRP78
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF: 1:100-500, P-ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 cell lysat...



Immunofluorescence analysis of COS7 cell...



Immunofluorescence analysis of Iovo cell...