

Product Datasheet

Acetyl-Histone H3(K27) antibody (orb2280251)

Description	Acetyl-Histone H3(K27) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthetic peptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Purity	≥95% as determined by SDS-PAGE
MW	15 kDa
Uniprot ID	P68431
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.