

## Product Datasheet

# AP conjugated Affinipure Rabbit Anti-Goat IgG (H+L) (orb2280267)

<b>Description</b>	AP conjugated Affinipure Rabbit Anti-Goat IgG (H+L)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	AP
<b>Tested Applications</b>	ELISA, WB
<b>Form/Appearance</b>	liquid
<b>Storage</b>	0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0. -20°C for 24 months , avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	1:1000-1:3000 for ELISA and Western blotting
<b>Expiration Date</b>	12 months from date of receipt.