

## Product Datasheet

# Sheep CXCL10/IP-10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit (orb2280278)

<b>Description</b>	Sheep CXCL10/IP-10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
<b>Reactivity</b>	Sheep
<b>Range</b>	31.2-2000pg/ml
<b>Target</b>	CXCL10/IP-10
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q95MZ7</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.