

## Product Datasheet

### Rat MT3 (Metallothionein 3) ELISA Kit (orb2280289)

<b>Description</b>	Rat MT3 (Metallothionein 3) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	MT3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.

