

## Product Datasheet

# Polyclonal Goat anti-rabbit IgG antibody (Goat anti-RB IgG) (orb2280479)

<b>Description</b>	Polyclonal Goat anti-rabbit IgG antibody (Goat anti-RB IgG)
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Purified rabbit IgG
<b>Preservatives</b>	50 mM NaCl 20 mM PBS pH 7.4 with 0.1% NaN <sub>3</sub>
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	2-8°C. For long-term storage, freeze at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Suitable for use in ELISA or Lateral flow immunoassay as a capture antibody. Other applications have not been tested but use in such assays should not be excluded. Great activity
<b>Clonality</b>	Polyclonal
<b>Purity</b>	> 95%
<b>Source</b>	Purified from whole goat anti-sera
<b>Expiration Date</b>	12 months from date of receipt.