

Product Datasheet

hnRNP K (Phospho-Ser284) antibody (orb228081)

Description Rabbit polyclonal antibody to hnRNP K

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Dilution Range WB: 1:500-1:2000

Expiration Date 12 months from date of receipt.