

---

## Product Datasheet

**Human CDCP1 Protein, Fc Tag (MALS verified)  
(orb2281239)**

**Description**

Human CDCP1 Protein, Fc Tag is

 (HEK293). It contains AA Phe 30 - Thr  
 ...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

CDCP1

**Preservatives**

PBS, pH 7.4

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 98.3 kDa. The protein migrates as 130-150 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**MW**

98.3 kDa

**Uniprot ID**
**Q9H5V8**
**Expiration Date**

6 months from date of receipt.