

Product Datasheet

Cynomolgus CCL1 Protein, Fc Tag (orb2281241)

Description	Cynomolgus CCL1 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Lys 24 - ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CCL1
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 34.9 kDa. The protein migrates as 40-43 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
MW	34.9 kDa
Uniprot ID	A0A8J8Y905
Expiration Date	6 months from date of receipt.