

---

## Product Datasheet

**Rat EphA2 Protein, His Tag (MALS verified)  
(orb2281243)**

**Description**

Rat EphA2 Protein, His Tag is

(HEK293). It contains AA Gln 24 - Met 5...

**Reactivity**

Rat

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

EphA2

**Preservatives**

PBS, pH 7.4

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.6 kDa. The protein migrates as 60-67 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**MW**

58.6 kDa

**Uniprot ID**
**D3ZBN3**
**Expiration Date**

6 months from date of receipt.