

Product Datasheet

Mouse Tenascin Protein, His Tag (orb2281247)

Description	Mouse Tenascin Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 23 - ...
Reactivity	Mouse
Conjugation	Unconjugated
Target	Tenascin
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 66.4 kDa.
MW	66.4 kDa
Uniprot ID	Q80YX1
Expiration Date	6 months from date of receipt.