

Product Datasheet

Human Tenascin-C Protein, His Tag (orb2281249)

Description	Human Tenascin-C Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 2 -...
Reactivity	Human
Conjugation	Unconjugated
Target	Tenascin-C
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.1 kDa.
MW	29.1 kDa
Uniprot ID	A0A994J5G0
Expiration Date	6 months from date of receipt.