

## Product Datasheet

# MOUSE ANTI HUMAN IgG (Fc) (CH2 DOMAIN):FITC ANTIBODY (orb2281262)

<b>Description</b>	MOUSE ANTI HUMAN IgG (Fc) (CH2 DOMAIN):FITC ANTIBODY
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Human IgG polyclonal.
<b>Preservatives</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Concentration</b>	IgG concentration 0.1 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Uniprot ID</b>	<b>P01861, P01859, P01860, P01857</b>
<b>Expiration Date</b>	12 months from date of receipt.