

Product Datasheet

MOUSE ANTI RABBIT CD45:FITC ANTIBODY (orb2281277)

Description	MOUSE ANTI RABBIT CD45:FITC ANTIBODY
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC
Immunogen	Glycoproteins isolated from the T cell line, RL-5.
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Expiration Date	12 months from date of receipt.