

## Product Datasheet

# MOUSE ANTI RABBIT CD8:RPE ANTIBODY (orb2281279)

<b>Description</b>	MOUSE ANTI RABBIT CD8:RPE ANTIBODY
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	RPE
<b>Tested Applications</b>	FC
<b>Preservatives</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
<b>Form/Appearance</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
<b>Storage</b>	+4°C DO NOT FREEZE
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Expiration Date</b>	12 months from date of receipt.