

## Product Datasheet

# Anti-Epstein-Barr Virus BFRF1 Antibody (orb2282879)

<b>Description</b>	Anti-BFRF1 Antibody catalog # orb2282879. This antibody reacts with Epstein-Barr Virus.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	500 µg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.